

# Sales Promotion

Valid till 31-03-2017








# 2 mag

magnetic *e* motion



Rs. 39,182

Order no: 30001  
**MIX 1**  
**Stirring volume: 1-10,000 L**  
 Speed range: 100 - 1,200 RPM  
 Stirring power (max): 20 W  
 Material housing: Stainless steel



Rs. 133,198

Order no: 30404  
**MIX 4 MS**  
**Stirring volume: 1 – 3,000 ml**  
 Speed range: 100 - 2,000 RPM  
 Stirring power: 28W (4x7)  
 Material housing: Stainless steel  
 Individual control over all stirring points




Rs. 19,500

Order no: 30101  
**MIX 1 eco**  
**Stirring volume: 1-3,000 L**  
 Speed range: 120 - 1,400 RPM  
 Stirring power (max): 10 W  
 Material housing: Stainless steel




Rs. 112,517

Order no: 30006  
**MIX 6**  
**Stirring volume: 1 – 3,000 ml**  
 Speed range: 100 - 2,000 RPM  
 Stirring power: 5/10/15/20 W  
 Material housing: Stainless steel



Rs. 205,482



Set up with a CD circulator, SS Bath tank and MIX6 magnetic stirrer combination

# BRANSON

# BREMA<sup>®</sup>

Ice Makers

2.8 L - Rs. 82,763  
5.7 L - Rs. 132,026



Order no: CPX-952-238R (2.8 liters)  
CPX-952-338R (5.7 liters)

### CPXH Series Ultrasonic baths

With sweep frequency  
40 kHz transducers  
99-minute digital timer  
Max heating temperature: 69 °C  
Perforated tray included!

Rs. 197,879



Order no: GB601A  
**GB 601 Ice flakes machine**  
Production in 24 hours: 60 kg  
Bin capacity: 10 kg  
Size: 450 x 620 x 680 mm (WXDXH)

# witeg

GERMANY

10 ml - Rs. 19,824  
50 ml - Rs. 28,974



Order no: 5370903 (10 ml)  
Order no: 5370905 (50 ml)  
**LABMAX Premium**  
Bottle top dispenser

10 ml: 0.20 division, ±0.06 tolerance and 0,010 deviation  
50 ml: 1.00 division, ±0.30 tolerance and 0,050 deviation



Rs. 57,139

Order no: 5497050  
**TITREX2000 Digital Burette**  
Volume: 0-50 ml  
Division: 0.01 ml  
Coeff. of Variation: 0.1%  
Accuracy: 0.2%

# CAT

Ingenieurbüro CAT  
M. Zipperer GmbH

Rs. 50,278

Order no: 60216-00  
**R 50D Overhead stirrer\***  
Volume: 100 L max  
Speed: 50 - 1600 RPM  
Torque: 55 Ncm  
Chuck: 0.5 to 10.5 mm  
Digital speed display  
\*Stirrer only



Order no: 60361-00  
**M21 Hotplate magnetic stirrer**  
Speed: 60-1600 RPM  
Round SS hot plate (Ø 135mm)  
Max 250 deg C  
Programmable switch off  
time delay



Rs. 54,561




Rs. 35,845

Order no: DH.WMH03021  
**MSH20D - Hotplate stirrer**  
Stirring volume: max 20 L  
Speed range: 80 - 1,500 RPM  
Temperature range: 380 °C, Heating power: 600 W  
Digital Feedback control,  
Ceramic coated Aluminium plate (180 x 180 mm)  
Digital LED display with Back-light function



Rs. 45,233

Order no: DH.WMH05020  
**MaXtir 500 - Hotplate stirrer**  
Stirring volume: max 20 L  
Speed range: 80 - 1,500 RPM  
Temperature range: 500 °C,  
Heating power: 850 W  
Digital Feedback control with Jog dial  
Solid ceramic glass plate (200 x 200 mm)  
Digital LED display with Back-light function



Rs. 78,890

Order no: DH.WMH03503  
**SMHS-3 - 3 point Hotplate stirrer**  
Stirring volume: max 20 L  
Speed range: 80 - 1,500 RPM  
Temperature range: 350 °C, Heating power: 400 W x 3  
Digital Feedback control,  
Ceramic coated Aluminium plate (140 x 140 mm)  
Digital LED display with Back-light function



Rs. 114,255

Order no: DH.WMH03506  
**SMHS-6 - 6 point Hotplate stirrer**  
Stirring volume: max 20 L  
Speed range: 80 - 1,500 RPM  
Temperature range: 350 °C, Heating power: 400 W x 6  
Digital Feedback control,  
Ceramic coated Aluminium plate (140 x 140 mm)  
Digital LED display with Back-light function



Rs. 40,966

Order no: DH.WHB00349  
**MaXtable H10 - Dry bath**  
Temperature range: Ambient + 5°C to 150°C, 300 W  
Timer: 99hr 59min, Digital LCD with Back-Light Function  
Heating block: 48 holes x for 1.5/2ml Centri. Tube

Metallic beads can also be supplied with MaXtable H10!  
Its another efficient way to use a dry bath without heating block!

Rs. 86,566





Rs. 10,882

Order no: DH.WCF00005


**CF-5 Mini-microcentrifuge**

RPM: 5,500 max, RCF: 1,775 x g

Capacity: 6x0.5-/1.5-/2.0-ml tubes,

2x0.2ml PCR 8-tube Strip or 16x0.2ml PCR 8 Tubes

Adapter: Included 6x0.2-/ 0.5-ml Tube Adapter



Rs. 88,865

Order no: DH.WSO04010

**SHO-1D Digital Orbital shaker**

Speed: 10-300 RPM, Max load: 8 kg

Timer: 99hr 59min (Continuous run)

Program: 10 Steps, 6 Programs

Universal Platform for various Shapes & Size of Flask

Capacity: 9 x 100 ml, 7 x 250 ml, 4 x 500/1000 ml Flask



Rs. 35,658

Order no: DH.WCF00010

**CF-10 Pro-microcentrifuge**

RPM: 1,000-13,500 max, RCF: 12,225 x g

Capacity: 0.2-/0.5-/1.5-/2.0-mlx12 Tubes

Rotor: Fixed Angle Rotor

Timer: 99hr 59min (Continuous run)



Rs. 13,656

Order no: DH.WVM00010

**VM-10 Vortex Mixer**

Speed: 0-3,300 RPM, Shaking: Orbital

Head: Pop-off cup and Ø76 mm platform

3-position Power Switch

Analogue phase controller



Rs. 297,538

Order no: DH.SWOV03030

**Thermostable SOV-30**

**Smart Vacuum Oven**

Range: 10~750 mm Hg, 0.63 mbar

Temperature: R/T +10 °C to 200 °C

Controller: Smart-Lab™ Controller

Shelves: 2

Internal : Stainless-steel ( # 304)

External : Powder Coated steel

Window : Tempered Safety Glass



Rs. 15,901



Order no: 1193  
**Dish shaped dewar flask**  
 With structured AL cover  
 Volume: 260 ml  
 Weight: 400 g

Rs. 18,954



Order no: 1194  
**Dish shaped dewar flask**  
 With structured AL cover  
 Volume: 400 ml  
 Weight: 500 g

Rs. 22,324



Order no: 1195  
**Dish shaped dewar flask**  
 With structured AL cover  
 Volume: 800 ml  
 Weight: 600 g



Rs. 13,268

Order no: 1068  
**Cylindrical Dewar Flask G9C type**  
 Blue coated metal cover incl. side grip  
 Volume: 10 L  
 Weight: 700 g  
 LN2 holding time: 60 h  
 Evaporation rate: 350 ml/day



Rs. 15,164

Order no: 1221  
**Dewar Carrying flask, Type: B-E**  
 With SS cover  
 Volume: 1 L  
 Weight: 1.5 Kg  
 LN2 holding time: 18 h  
 Evaporation rate: 1,2 L/day

**Terms and Conditions:**

Price Basis:	Ex-works, Bangalore
Sales Tax:	VAT/CST @ 14.5% extra
Payment:	100% Advance by RTGS/NEFT
Bank Details:	Beneficiary - Scientific Research Instruments Co. Put, Ltd., Bank name: ICICI Bank, M. G. Road Branch, Bangalore A/c. No. 000205023122 SWIFT CODE: ICICI NBBCTS NEFT CODE: ICIC0000002
Delivery:	2-4 Weeks from the date of Purchase Order
Delivery Terms:	F.O.R Destination for End Users F.O.R Bangalore for Distributors
Warranty:	1 year from the date of invoice
Validity:	Until stock lasts and please check availability of stock with head office before committing to customer